# FITAGE

## USER MANUAL

MANUAL LANGUAGES:

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**SCAN THE QR CODE**

to view the digital manual  
in all languages

# WELCOME TO FITAGE GET IN SHAPE

Welcome to FITAGE, where we believe that everyone has the power to take control of their health and wellbeing. We are thrilled to introduce you to our advanced blood pressure monitor, which will be your personal health assistant on your journey towards better health.

With clinically approved measuring techniques and an easy-to-use design, our monitor provides accurate measurements that you can rely on. Our custom-made large LED-screen ensures that you can read your results quickly and clearly, making tracking your health a breeze.

At FITAGE, we understand that monitoring your health can be daunting, but we believe that every small step you take towards better health is worth celebrating.

Whether you are looking to manage an existing condition or simply want to stay on top of your health, our monitor is here to support you every step of the way.

So, are you ready to take the first step towards a healthier you? Let FITAGE be your guide and start your journey today.

**SAFETY NOTICE.** Intended Use: The digital blood pressure monitor is designed to measure the systolic and diastolic blood pressure as well as the pulse rate of adults and adolescents over the age of 12. This is done using a non-invasive technique on the upper arm. The device is suitable for both clinical and home use, and all functions can be safely utilized. Measurement values can be easily read on the LCD display. It's important to note that the device should only be used on the upper arm of adult patients. The individual operating the device is considered the intended operator.

# SPECIFICATIONS

<b>Product Description</b>	FITAGE Arm-type Blood Pressure Monitor	
<b>Model</b>	FITAGE-06	
<b>Display</b>	LCD Digital Display Oscillometric Method	
<b>Measurement Method</b>	Size: 102.1 mm × 68.9 mm (4.02" x 2.71")	
<b>Measurement Range</b>	Systolic Pressure	60mmHg~260mmHg
	Diastolic Pressure	30mmHg~200mmHg
	Pressure range	0mmHg~300mmHg
	Pressure deviation	±3mmHg
	Pulse	30 ~ 180 Beats/Minute
	Pulse deviation	±5%
<b>Pressurization</b>	Automatic Pressurization	
<b>Memory</b>	120 Memories in Two Groups with Date and Time	
<b>Function</b>	Irregular Heartbeat Detection	
	WHO Classification Indicator	
	Test Average over Three Measurements	
	Low Battery Detection	
	Automatic Power-Off	
	Voice	
	Backlight	
<b>Power Source</b>	4x AA batteries and USB-C cable	
<b>Battery Life</b>	Approximately 2 months at 3 tests per day	
<b>Unit Weight</b>	Approximately 324 g (11.43 oz.) (excluding batteries)	
<b>Unit Dimensions</b>	Approximately 166 x 114 x 72mm (6.53" x 4.49" x 2.83" ) (L x W x H)	
<b>Cuff Circumference</b>	Medium/large cuff: Fits arm circumference 22-42 cm	

<b>Operating Environment</b>	Temperature	10°C ~ 40°C (50°F~104°F)
	Humidity	15% ~ 93%RH
	Pressure	700hPa~1060hPa
<b>Storage Environment</b>	Temperature	-25°C~70°C (-13°F~158°F)
	Humidity	≤93% RH
<b>Classification:</b>	Internal Powered Equipment, Type BF  , Cuff is the Applied Part	
<b>Ingress Protection Rating:</b>	IP20, Indoor Use Only	

The specifications provided may change at any time without prior notice.

It's good to know that the Blood Pressure Monitor meets European regulations and has passed various safety and performance standards. This ensures that the device is safe and reliable for use in measuring blood pressure.

The Blood Pressure Monitor bears the CE mark "CE 0123" and complies with European regulations. Additionally, it adheres to various safety, EMC, and performance standards, including but not limited to: EN 60601-1 for general safety requirements, EN 60601-1-2 for electromagnetic disturbance requirements, IEC80601-2-30 for essential performance and safety of automated non-invasive sphygmomanometers, EN 1060-3 for supplementary requirements of electro-mechanical blood pressure measuring systems, and ISO 81060-2 for clinical validation of automated measurement type of non-invasive sphygmomanometers.

## Correct Disposal of This Product

(Waste Electrical & Electronic Equipment)

The product contains a marking indicating that it should not be disposed of with regular household waste when it reaches the end of its useful life. To avoid potential harm to the environment or human health, it is important to separate this product from other types of

waste and recycle it responsibly. To dispose of this product in an environmentally safe recycling center, contact the retailer where the product was purchased or your local government office for instructions. Business users should contact their supplier and review the purchasing agreement terms and conditions. Do not mix this product with other commercial waste for disposal. This product does not contain hazardous materials.

## ELECTROMAGNETIC COMPATIBILITY INFORMATION

This device meets the electromagnetic compatibility (EMC) requirements outlined in the international standard IEC 60601-1-2. However, it is important to follow the specific conditions listed in the table below in order to satisfy these requirements. Since the device is an electrical medical product, it is subject to special precautionary measures with regard to EMC which will be provided in the instructions for use. Please note that the device may be affected by portable or mobile HF communications equipment, and using non-approved accessories with it can have a negative impact on its EMC. Additionally, it is recommended to avoid using the device in close proximity to other electrical equipment.

**Table 1**

The manufacturer provides guidance on the electromagnetic emissions of the device.		
This device is intended to be used only in a specific electromagnetic environment, which is specified below. It is the responsibility of the customer or user to ensure that the device is used only in this specified environment.		
Emissions test	Compliance	Electromagnetic environment -guidance
<b>Radiated emission CISPR 11</b>	Group 1, class B.	This device only uses RF energy for its internal operations, which means that its emissions are minimal and unlikely to interfere with other electronic equipment nearby.

<b>Conducted emission CISPR 11</b>	Group 1, class B.	This device can be used in any type of establishment, including residential properties and buildings that are directly connected to the public low-voltage power supply network that are used for residential purposes.
<b>Harmonic emissions IEC 61000-3-2</b>	Class A	
<b>Voltage fluctuations/ flicker emissions IEC 61000-3-3</b>	Complies	

**Table 2**

The manufacturer provides guidance on the electromagnetic immunity of the device.			
This device is designed to be used in a specific electromagnetic environment, which is specified below. It is the responsibility of the customer or user to ensure that the device is used only in this designated environment.			
<b>IMMUNITY test</b>	<b>IEC 60601 test level</b>	<b>Compliance level</b>	<b>Electromagnetic environment -guidance</b>
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	± 8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	Floors should be made of wood, concrete, or ceramic tile. If the floors are covered with synthetic material, the relative humidity must be maintained at a minimum of 30%.
Electrostatic transient/ burst IEC 61000-4-4	± 2 kV , 100kHz, for AC power port	± 2 kV , 100kHz, for AC power port	The quality of the main power supply should be consistent with what is typically found in a commercial or hospital setting.
Surge IEC 61000-4-5	±0.5kV, ±1kV (differential mode)	±0.5kV, ±1kV (differential mode)	The quality of the main power supply should be consistent with what is typically found in a commercial or hospital setting.

Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	% UT; 0,5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°  0 % UT; 1 cycle and 70 % UT; 25/30 cycles Single phase: at 0°, 0% UT; 250/300 cycle	% UT; 0,5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°  0 % UT; 1 cycle and 70 % UT; 25/30 cycles Single phase: at 0°, 0 % UT; 250/300 cycle	The quality of the main power supply should be consistent with what is typically found in a commercial or hospital setting.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m; 50Hz or 60Hz	30 A/m; 50Hz or 60Hz	The levels of power frequency magnetic fields should match what is typically found in a commercial or hospital setting.

**Table 3**

The manufacturer provides guidance on the electromagnetic immunity of the device.			
This device is designed to operate within a specific electromagnetic environment, which is specified below. It is the responsibility of the customer or user to ensure that the device is used only in this designated environment.			
IMMUNITY test	IEC 60601 test level	Compliance level	Electromagnetic environment -guidance
Conducted RF IEC 61000-4-6	3V for 0.15- 80MHz; 6V in ISM and amateur radio bands between 0.15- 80MHz	3V for 0.15- 80MHz; 6V in ISM and amateur radio bands between 0.15- 80MHz	When using portable or mobile RF communication equipment, you should maintain a certain distance from the device, including its cables. This distance should be calculated based on an equation that corresponds to the frequency of the transmitter, and it's recommended that you do not use the equipment closer than this distance. Recommended separation distance
Radiated RF IEC 61000-4-3	385MHz, 27V /m	385MHz, 27V /m	

			$d = \left[ \frac{3.5}{E_1} \right] \sqrt{P}$ $d = \left[ \frac{7}{E_1} \right] \sqrt{P}$ <p>80 MHz to 800 MHz</p> <p>800 MHz to 2.7 GHz</p> <p>The equation for calculating the recommended separation distance between you and the transmitter depends on two factors: the maximum output power rating of the transmitter in watts (W), denoted by "P," and the distance in meters (m), denoted by "d."</p> <p>The strength of electromagnetic fields emitted by fixed RF transmitters should be measured through an electromagnetic site survey. The measurements must show that the field strengths are below the compliance level for each frequency range.</p> <p>If you see this symbol on equipment, there is a possibility that interference may occur in its vicinity:</p> 
	450MHz, 28V /m	450MHz, 28V /m	
	710MHz,745 MHz,780MHz 9V/m	710MHz,745 MHz,780MHz 9V/m	
	810MHz,870 MHz,930MHz 28V/m	810MHz,870 MHz,930MHz 28V/m	
	1720MHz,1845 MHz, 1970MHz 28V/m	1720MHz,1845 MHz, 1970MHz 28V/m	
	5240MHz,5500 MHz,5785MHz 9V/m	5240MHz,5500 MHz,5785MHz 9V/m	

Table 4

The following are the recommended distances that should be maintained between portable and mobile RF communications equipment and the device:		
The device is designed to be used in an electromagnetic environment that controls radiated disturbances. To avoid electromagnetic interference, users should maintain a minimum distance between the device and portable or mobile RF communication equipment (transmitters). The recommended distance depends on the maximum output power of the communication equipment.		
Rated maximum output power of transmitter.	The recommended separation distance varies based on the frequency of the transmitter and is measured in meters (m).	
	80 MHz to 800 MHz $d = \left[ \frac{3.5}{E_1} \right] \sqrt{P}$	800 MHz to 2.7 GHz $d = \left[ \frac{7}{E_1} \right] \sqrt{P}$
W		
0.01	0.12	0.23
0.1	0.38	0.73
1	1.2	2.3
10	3.8	7.3
100	12	23
<p>If the maximum output power of a transmitter is not listed in the previous guidelines, you can estimate the recommended separation distance (d) in meters (m) using an equation that corresponds to the frequency of the transmitter. To use the equation, you need to know the maximum output power rating of the transmitter in watts (W), as specified by the manufacturer.</p> <p>Note 1: When it comes to separation distance, if the frequency of the transmitter is either 80 MHz or 800 MHz, you should follow the guidelines for the higher frequency range.</p> <p>Note 2: Please note that these guidelines may not be applicable in all circumstances since electromagnetic propagation can be influenced by factors such as absorption and reflection from structures, objects, and people.</p>		

# ADDITIONAL NOTES

Important Instructions Before Use:

**Warning:** Avoid using this equipment near or stacked with other devices to prevent malfunction. If necessary, ensure both this equipment and others operate normally when used together.

**Warning:** Keep portable RF communication equipment, including antennas and cables, at least 30 cm (12 inches) away from any part of the Fully Automatic Digital Blood Pressure Monitor to avoid performance issues.

The software identifier is referenced in the software evaluation report, with the file code NYRJ211221001.

**Verify Manometer Pressure Accuracy:** Follow these steps to check the accuracy of the pressure:

While the device is off, press and hold the "START/STOP" button and insert batteries.

Release the "START/STOP" button when the LCD screen is full.

When "00" appears on the screen, the device is static.

Connect a calibrated standard pressure gauge and manual pressure device to the monitor.

Apply manual pressure within the display range and compare readings with the standard gauge.

**Contraindications:** This product is not suitable for infants or individuals who cannot communicate.

**Patient as Operator:** Patients should only replace batteries and not perform other maintenance tasks.

**Warning:** Unauthorized modification of this equipment is prohibited.

**Essential Performance Maintenance:** Calibration is performed before leaving the factory. If accuracy deviates significantly, contact the manufacturer for recalibration.

**Mechanical Strength and Heat Resistance:** The device maintains heat resistance throughout its service life.

**Warning:** Improper placement of the monitor and cuff can be hazardous, especially to children.

**Biocompatibility:** Both cuff and monitor case are tested for biocompatibility and should be discontinued if allergic reactions occur.

**Warning:** Only professionals should modify the equipment to ensure accurate measurements.

**Warning:** Prolonged exposure to harsh conditions may impair equipment performance.

**Warning:** This device is not suitable for use on children or pets.

**Cleaning:** Lay operators should follow the provided cleaning instructions.

**Warning:** Do not use a damaged cuff for blood pressure measurement.

**Warning:** If discomfort occurs during measurement, deflate the cuff or remove it promptly.

**Warning:** If unexpected readings occur, take multiple measurements and consult a doctor.

**Warning:** Using the equipment outside the specified environment may damage it and lead to inaccurate readings.

**Warning:** Loss or degradation of essential performance due to electromagnetic disturbances should be reported to customer service.

**Warning:** Use only specified accessories and cables to avoid equipment malfunction.

**Warning:** Failure to use the equipment in a shielded location may degrade performance or interfere with other devices.

**Warning:** This equipment is not suitable for use with flammable agents or in oxygen-rich environments.

**Indications for Use:**

The digital blood pressure monitor measures blood pressure and pulse rate from the upper arm.

### **Correct Disposal of This Product:**

Dispose of this product responsibly by separating it from other wastes and recycling it. Contact your retailer or local government office for disposal instructions. Business users should refer to their purchasing agreement for disposal guidelines. This product does not contain hazardous materials.


# AFTER-SALE SERVICE

After purchasing our product, if you need any assistance, please contact our after-sale service, FITAGE B.V., located at Kabelweg 80, 2371DX Roelofarendsveen, The Netherlands. You can also reach out to them via email at **[support@fitage.nl](mailto:support@fitage.nl)** or visit their website at **[www.fitage.nl](http://www.fitage.nl)** for any support or assistance. Our dedicated team at FITAGE B.V. will be more than happy to assist you with any questions or concerns you may have about our product.

# WARRANTY



The Blood Pressure Monitor comes with a 2-year warranty from the date of purchase. In case of any malfunction due to defective components or poor workmanship, we will repair or replace the unit free of charge. However, please note that the warranty does not cover damages caused by improper handling of the product. For more information, please contact your local retailer or our warranty policy on **[www.fitage.nl](http://www.fitage.nl)**.

# DEVICE ILLUSTRATION


# DEVICE ILLUSTRATION

Arm Cuff

# DEVICE ILLUSTRATION


## Contents


1. Monitor

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

3. Arm Cuff

4. Storage Case

ZMOPPMedical AC Adapter  
(recommended, not provided)

5. USB-C Cable

# IMPORTANT GUIDELINES FOR BLOOD PRESSURE TESTING

1. Refrain from eating, exercising, and taking a bath for at least 30 minutes before testing your blood pressure.
2. Sit in a peaceful environment for a minimum of 5 minutes before testing.
3. Avoid standing while testing. Sit comfortably with your arm at heart level.
4. Refrain from talking or moving while testing.
5. Avoid strong electromagnetic interference, such as microwave ovens and cell phones, while testing.
6. Wait for at least 3 minutes before re-testing.
7. Try to measure your blood pressure at the same time every day to maintain consistency.
8. Only compare test results when the monitor is used on the same arm, in the same position, and at the same time of day.
9. This blood pressure monitor is not recommended for individuals with severe arrhythmia.
10. Do not use this blood pressure monitor if the device is damaged.

# QUICK START GUIDE

To get accurate readings, avoid smoking, eating, drinking caffeinated drinks, or exercising for 30 minutes before taking the measurement.

1. Install batteries as shown in Figure A.
2. Insert the cuff air plug into the left side of the monitor as shown in Figure B.


Figure A


Figure B

3. Remove any thick clothing from your arm.
4. Rest for several minutes before testing. Sit down in a quiet place, preferably at a desk or table, with your arm resting on a firm surface. Keep your legs uncrossed and your feet flat on the floor, as shown in Figure C.




Figure C

5. Apply the cuff to your left arm and make sure it is level with your heart. The bottom of the cuff should be placed approximately 1-2 cm (0.4-0.8") above your elbow joint, as shown in Figures D and E.


6. Press "Start/Stop" button to start testing.

## OPERATION

Instructions for Installing Batteries:

1. Slide the battery cover off in the direction of the arrow.
2. Install four new AA alkaline batteries according to the polarity (+/-) shown inside the battery compartment.
3. Close the battery cover.

You can find the USB-C adapter jack on the back of the monitor. It's recommended to the provided USB-C cable. When connecting the adapter, make sure that the pin is positive inside and negative outside. Do not use any other type of adapter as it may damage the device.

USB-C Adapter Jack

Note: The power supply is specified as part of ME equipment.


# SYSTEM SETTINGS

System Settings:

To access the system settings, turn off the device if needed (by pressing Start/Stop) and press the "SET" button. The Memory Group icon will start flashing.

## 1. Select Memory Group

While in the System Setting mode, you can store test results into two different groups. This allows multiple users to save individual test results (up to 60 memories per group). To choose a group setting, press the "M" button. The test results will automatically save in the selected group.

## 2. Setting Time and Date

To set the time and date, press the "SET" button a couple of times while the device is powered off until the Time/Date mode is activated. Use the "M" button to set the year, then press the "SET" button to confirm the current year. Continue setting the date, hour, and minute by using the "M" button and confirming each selection by pressing the "SET" button. Each time the "SET" button is pressed, it will save your selection and move on to the next setting (month, day, hour, and minute).

### 3. Time Format Setting

To adjust the time format, press the "SET" button once more to enter the Time Format setting mode. Then, use the "M" button to choose between the two options: EU for European Time or US for U.S Time. Once you've made your selection, press the "SET" button again to confirm and lock it in.

### 4. Voice Setting

To adjust the voice setting, press the "SET" button to enter the voice setting mode. Then, use the "M" button to turn the voice format ON or OFF as desired.

## 5. Volume Setting

To adjust the volume of the voice feature, press the "SET" button to enter the volume setting mode. Then use the "M" button to select from six different volume levels.

## 6. Saved Settings

You can save any changes you make to the system settings by turning off the device. To turn off the device, press the "START/STOP" button while in any setting mode. The device will automatically save all information if it is left on and not in use for 3 minutes.

## 7. Backlight Setting

Press "SET" button to enter backlight setting mode. Set backlight format ON or OFF by pressing the "M" button.


## 8. Audio Function


- Press "SET" to enter the menu
- Keep pressing "SET" to navigate to the language settings (L)
- Press "MODE" to select a language:
  - English: L1
  - Dutch: L2
  - German: L3

# APPLY THE ARM CUFF

1. Insert the air plug firmly into the opening on the left side of the monitor.

2. Make sure the sticky nylon section is facing outward, then insert the end of the cuff underneath the metal ring of the cuff.
3. Position the cuff about 1-2 cm (0.4-0.8 inches) above the elbow joint and fasten it securely. For optimal results, apply the cuff to a bare arm and keep it level with your heart while performing the test.


Please do not insert the air plug into the opening on the right side of the monitor as it is only intended for an optional power supply.


## MEASUREMENT

### 1. Turn the Device On

To turn on the device, press and hold the "START/STOP" button. The LCD screen will appear for one second while the device performs a quick diagnosis. Once it is ready for testing, a voice tone will indicate that it is ready.

It's important to note that the device will not work if there's still air from the previous test in the cuff. If this is the case, the LCD screen will flash ↓ until the pressure is stable.


## 2. Adding Pressure

The device will automatically inflate the cuff to the correct pressure and stop inflating. It's important to remain quiet during this process.

Note: If the cuff is not properly applied to the arm, pressurization will gradually subside and ultimately stop. In this case, press the "START/STOP" button to turn off the device.


## 3. Testing

After the cuff inflates, the air pressure will slowly decrease, and the corresponding cuff pressure value will be displayed on the screen. At the same time, a flashing symbol will appear on the screen indicating that the heart-beat is being detected.

Note: It is important to remain relaxed during the testing process. Avoid speaking or moving any part of your body.

## 4. Result Display

After the testing is complete, the device will display the measurements for systolic and diastolic blood pressure, and a voice broadcast will announce the results. An indicator representing the current measurement will appear next to the corresponding WHO classification.

Note: Please refer to pages 23-24 for detailed information on WHO Blood Pressure Classification.

### **Irregular Heartbeat Indicator**

The Irregular Heartbeat Indicator is a feature that alerts you if the monitor detects an irregular heart rhythm two or more times during the measuring process. When this happens, the Irregular Heartbeat Symbol “♥” will appear on the screen along with your measurement results. An irregular heartbeat rhythm is defined as a rhythm that is either 25% slower or faster than the average rhythm detected while measuring systolic and diastolic blood pressure. If you frequently see the Irregular Heartbeat Symbol “♥” with your test results, it is recommended that you consult your physician.

### **5. Deleting and Storing Test Results**

If you are not satisfied with your test result or have any other reason to delete it, you can do so by pressing the “SET” button right after the result is displayed. If you don't delete it, it will be automatically saved by date in the Memory Group you have selected before testing. Please make sure to choose the correct Memory Group before testing.

The device can store up to 120 test results per Memory Group. If the number of tests exceeds 120, the oldest readings will be deleted to make room for the new ones, and the most recent results will be displayed first.


# POWER OFF

You can turn off the device in any mode by pressing the "START/STOP" button. The device will turn off itself automatically after 3 minutes of no operation in any mode.


**Safety Precaution:** If the pressure in the arm cuff becomes too extreme during testing, press the "START/STOP" button to turn off the power. The cuff pressure will rapidly decrease once the device is turned off.

# MEMORY CHECK AND TEST AVERAGE

When the power is off, you can check your past test results and the average using the "M" button. Pressing the "M" button for the first time will display the average value of all saved measurements from the selected group with a voice broadcast (shown as AVGALL). Pressing the "M" button again will show the average of the last seven days' test results between 5:00 and 8:59 in the morning with a voice broadcast (shown as AVGM). Pressing the "M" button again will show the average of the last seven days' test results between 18:00 and 19:59 in the evening with a voice broadcast (shown as AVGP). Note that all times are in European format.


To view previous test results, press the "M" button again. You can view the most recent test result. When the test results are displayed, you can press the "M" button to browse through all the test results stored in memory. To view previous tests in other memory groups, you must first select the desired group and then turn off the monitor (refer to "Select Memory Group" on page 10).

Note: If there are no saved results, the systolic blood pressure, diastolic blood pressure, and pulse rate will display as "...".


## DELETE MEASUREMENT

To delete the measurement records from a selected group, go to Memory Check mode and hold the "SET" button for about 3 seconds. You will hear a voice broadcast saying "Memory Clear", and then the device will switch to testing mode. After this, press the "START/STOP" button to turn off the device.

Note: Once record data has been deleted, it cannot be recovered.

# LOW BATTERY INDICATOR


When the battery life is low and unable to inflate the cuff for testing, the device will announce "Low Battery". At the same time, the "🔋" symbol will appear on the screen for about 5 seconds before turning off. Replace the batteries at this time. This process will not result in any memory loss.

# STATIC PRESSURE MEASUREMENT

To measure static pressure, follow these steps:

1. Turn off the blood pressure meter.
2. Press and hold the "START/STOP" button and insert the batteries.
3. Keep holding the "START/STOP" button until the LCD screen is completely filled.
4. Release the "START/STOP" button when the screen shows double zero.
5. The blood pressure meter is now in static state, and you can see the software version on the heart rate display.

Note:

Only service personnel are permitted to access this mode.

This mode is unavailable for normal use.

# TROUBLESHOOTING

Problem	Possible Cause	Solution
The blood pressure results are not within a typical range	The cuff is too tight or not properly positioned on the arm	Reposition the cuff approximately 1-2 cm (1/2") above the elbow joint
	Inaccurate test results may be due to body movement or monitor movement during testing.	To ensure accurate results, sit in a relaxed position with your arm placed near your heart. Avoid speaking or moving your body parts while testing. Additionally, make sure the monitor device is placed in a stationary position throughout the testing period.
"ERR" is displayed	The cuff fails to inflate properly	Make sure the hose is properly fastened to the cuff and monitor device
	Improper operation	Read the user manual carefully and re-test properly.
	Pressurization is over cuff rated pressure 300mmHg	Read the user manual carefully and re-test properly.
	Values of blood pressure that fall outside the specified range for blood pressure.	Please re-test properly.
"  "	Appear when the batteries life is depleting and unable to inflate cuff for testing	Replace the batteries as soon as the symbol appears.

## BLOOD PRESSURE INFORMATION


### Blood Pressure

Blood pressure refers to the pressure exerted by blood against the walls of arteries, usually measured in millimeters of mercury (mmHg). The highest pressure exerted against blood vessel walls when the heart beats is called systolic blood pressure, while the pressure exerted on blood vessels when the heart is resting between beats is called diastolic blood pressure.

Blood pressure can vary throughout the day due to different factors. Excitement, tension, alcohol consumption, and bathing can all affect blood pressure. Hormones such as adrenaline, which is

released during stress, can cause blood vessels to narrow, leading to an increase in blood pressure.


When blood pressure readings are too high, it means the heart is working harder than it should.


Example: fluctuation within one day (male, 35 years old)

## WHO BLOOD PRESSURE CLASSIFICATION INDICATOR

The blood pressure monitor has a feature called the WHO Blood Pressure Classification Indicator which follows the guidelines from the World Health Organization. The monitor is color-coded to show the test results and their classification based on the chart provided.


Blood Pressure Classification Indicator

## HEALTH REMINDER

It's important to take hypertension seriously as it can have serious consequences on your health and quality of life. It can lead to issues such as heart failure, kidney failure, and cerebral hemorrhaging. However, by adopting a healthy lifestyle and regularly visiting your doctor, hypertension and related diseases can be managed more easily if diagnosed early on. Remember to prioritize your health and well-being by taking the necessary steps to prevent and manage hypertension.


It's important not to panic if you get an unusual blood pressure reading. It's recommended to take 2-3 readings at the same time each day over a period of time to get a better idea of your blood pressure. If your readings continue to show abnormalities, it's advised to consult with your physician.

## MORNING HYPERTENSION

Morning hypertension refers to high blood pressure that occurs specifically in the morning hours, typically within one to two hours of waking up. This is different from the normal fluctuations in blood pressure that occur throughout the day.

Morning hypertension is a risk because it has been associated with an increased risk of cardiovascular events, such as heart attacks and strokes. This is thought to be due to the fact that blood pressure is naturally higher in the morning, and this additional increase can put extra strain on the heart and blood vessels. Additionally, morning hypertension can lead to an increased risk of other conditions such as kidney disease and cognitive decline. It is important for people with hypertension to monitor their blood pressure regularly and discuss any concerns with their healthcare provider.

The Symbol for Morning Hypertension is when on the LCD the up arrow and the time both flash at the same time. This symbol is displayed if the morning average reading for a week is above 140 for the Systolic Blood Pressure value and/or 90 for the Diastolic Blood Pressure value when viewing the Morning Averages (AVGAM) and the Afternoon Averages (AVGPM) in the memory. You can view measurement averages by pressing the (M) button (see Memory Check and Test Averages).

# BLOOD PRESSURE Q&A

## **Q: How does measuring blood pressure at home differ from measuring it at a professional healthcare clinic?**

A: Measuring blood pressure at home is now considered to be more accurate as it reflects an individual's daily life. Blood pressure readings may be elevated in a clinical setting due to anxiety or nervousness, known as White Coat Hypertension.

Note: Abnormal test results may be caused by:

1. Improper cuff placement - ensure the cuff is snug, not too tight or too loose, and positioned approximately 1-2 cm (1/2") above the elbow joint.
2. Improper body position - keep your body in an upright position.
3. Anxiety or nervousness - take 2-3 deep breaths, wait a few minutes, and resume testing.

## **Q: What are the factors that can cause different blood pressure readings?**

A: Blood pressure can vary throughout the day due to factors such as diet, stress, cuff placement et cetera.

## **Q: Should I use my left or right arm for measuring blood pressure?**

A: Either arm can be used for testing, but for accurate comparisons, the same arm should be used consistently. Testing on the left arm may provide more accurate results as it is closer to the heart.

## **Q: What is the best time of day to measure blood pressure?**

A: The morning or any time when you are relaxed and stress-free is ideal for testing.

# MAINTENANCE


1. Do not drop, slam, or throw the device to prevent damage.

2. Keep the device away from extreme temperatures and direct sunlight.

3. Use a soft cloth and mild detergent to clean the device. Wipe off any excess detergent with a damp cloth.


4. To clean the cuff, use a soft cloth with rubbing alcohol and then wipe it with a damp cloth. Let the cuff dry at room temperature. Remember to clean and disinfect the cuff between different users.

5. Do not use petrol, thinners, or similar solvents.

6. Remove batteries when the device is not in use for an extended period of time.

7. Do not attempt to disassemble the product.

8. It is recommended to have the device's performance checked every two years.


9. The expected service life of the device is approximately three years with 10 tests per day.

10. Only service personnel should perform maintenance on the device. If parts or technical support are needed, they will be provided.

# PRECAUTIONS AND SAFETY MEASURES

This device is designed for adults and adolescents aged 12 years or older to measure their systolic and diastolic blood pressure and pulse rate. All of the device's functions can be used safely, and the readings can be viewed on the LCD display. The device should only be placed on the upper arm of an adult when taking a measurement. Before using the device, it's essential to read the manual thoroughly and keep it for future reference. If you have any questions or concerns about your blood pressure, consult your doctor.

To prevent any risks or damage, please follow all of the warning precautions provided in the manual. Use the device only as intended and read all of the instructions before using it.

WARNING SIGNS AND SYMBOLS USED	
	Caution
	Mandatory
	Prohibited
	Type BF (body floating) equipment
	Instructions before use must be consulted
	Serial Number
	Discard the used product to the recycling collection point according to local regulations
	The product conforms to the requirements of the Regulation (EU) 2017/745 MDR on medical devices
	Manufacturer
	Authorized representative within the European Union
	Keep dry
	Keep out of sunlight
	Manufacturing date
	Medical Device

### CAUTION

If you have serious circulation problems, you may experience discomfort while using this product. It's important to consult your physician before using it.

If your test results consistently show abnormal readings, contact your physician. Do not attempt to self-treat these symptoms without first consulting with your physician.

This product is only intended to be used for its designed purpose. Do not misuse it in any way.

This product is not meant for infants or individuals who are unable to communicate their intentions.

Do not attempt to disassemble or repair the product.

Do not use cell phones or other devices that generate strong electrical or electromagnetic fields near the device, as this may cause incorrect readings and interfere with the device.

Use only a recommended USB-C adapter that complies with EN 60601-1 and EN 60601-1-2 (see page 6). Using an unauthorized adapter may cause fire or electric shock.

### USB-cable advice

Before using the device, please read the user's manual carefully and use the USB cable that is specified in the manual.

Do not insert or remove the USB cable while the device is transferring data. To replace the batteries of the device, first remove the USB cable.

Do not touch, insert, or remove the USB cable with wet hands. When the device is not in use, disconnect the USB cable from the device and do not leave it connected.

Battery Precautions
Do not use new and old batteries together in the device at the same time.
Replace batteries when low battery indicator  appears on the screen.
Make sure to insert the batteries correctly with the correct polarity.
Use only one type of battery in the device. Long-life alkaline batteries are recommended.
If the device will not be used for more than three months, remove the batteries.
Dispose of the batteries in a proper manner in accordance with local laws and regulations.

## SAFETY NOTICE

Important Instructions before Use:

1. Self-monitoring is not the same as self-diagnosis. A healthcare professional who knows your medical history should interpret your blood pressure measurements.
2. If your test results show abnormal readings consistently, contact your physician.
3. If you take medication, consult your physician to determine the best time to measure your blood pressure. Do not change your prescribed medication without consulting your physician.
4. If you have circulation problems, consult your physician before using the device.
5. For people with irregular or unstable circulation caused by diabetes, liver disease, arteriosclerosis, or other medical conditions, blood pressure values measured at the wrist may differ from those measured at the upper arm. Nevertheless, monitoring the trends in your blood pressure taken at either the arm or the wrist is useful and important.
6. People suffering from vascular constriction, liver disorders or diabetes, people with cardiac pacemakers or a weak pulse, and pregnant women should consult their physician before measuring their blood pressure themselves. Their condition may produce different values.

7. People with arrhythmias such as atrial or ventricular premature beats or atrial fibrillation should only use this blood pressure monitor in consultation with their doctor. In certain cases, the oscillometer measurement method can produce incorrect readings.
8. Frequent measurements can cause injury to the patient due to blood flow interference.
9. Do not apply the cuff over a wound as this can cause further injury.
10. Do not attach the cuff to a limb being used for IV infusions, any other intravascular access, therapy, or an arterio-venous (A-V) shunt. The cuff inflation can temporarily block blood flow, potentially causing harm to the patient.
11. Do not place the cuff on the arm on the side of a mastectomy. In the case of a double mastectomy, use the side of the least dominant arm.
12. Pressurization of the cuff can temporarily cause loss of function of simultaneously used monitoring equipment on the same limb.
13. A compressed or kinked connection hose may cause continuous cuff pressure resulting in blood flow interference and potentially harmful injury to the patient.
14. Make sure the device does not result in prolonged impairment of the circulation of the patient.
15. Only use the product for its intended purpose. Do not misuse it in any way.
16. The product is not intended for infants or individuals who cannot express their intentions.
17. Prolonged over-inflation of the bladder may cause ecchymoma of the arm.
18. Do not disassemble or attempt to repair the device or arm cuff.
19. Only use the approved arm cuff for this device. Using other arm cuffs may result in incorrect measurement results.
20. The system may produce incorrect readings if stored or used outside the manufacturer's specified temperature and humidity ranges. Make sure to store the blood pressure monitor in a place where children, pets, and pests cannot access it.

21. Do not use the device near strong electrical or electromagnetic fields generated by cell phones or other devices, as they may cause incorrect readings and interference, or become an interference source to the device.
22. Do not mix new and old batteries simultaneously.
23. Replace both batteries when Low Battery Indicator  appears on screen.
24. Use only long-life alkaline batteries for this device.
25. Remove the batteries from the device when not in use for more than 3 months.
26. Insert the batteries with their polarities correctly aligned.
27. Dispose of batteries properly and observe local laws and regulations.
28. Only use the recommended USB-C cable that complies with EN 60601-1 and EN 60601-1-2. An unauthorized adapter may cause fire and electric shock.
29.  Please refer to the instruction manual/booklet for guidance on operating the device.
30. The PC that connects to the device via USB must meet the requirements of either IEC 60601-1 or IEC 60950-1 standard.
31. To ensure measurement accuracy, do not use the device during transport vehicles such as ambulances or helicopters.
32. This device contains small parts that may be a choking hazard for infants if swallowed.
33. When replacing the batteries, please align the polarities of each battery with the +ve and -ve signs imprinted on the battery housing.
34. The time required for the device to warm up from the minimum storage temperature (-25°C) until it's ready for use at ambient temperature (20°C) is about 2 hours.
35. The time required for the device to cool down from the maximum storage temperature (70°C) until it's ready for use at ambient temperature (20°C) is about 2 hours.

## WE WANT TO LISTEN TO OUR CUSTOMERS

Your opinion is everything for FITAGE. Your experience drives us to help contributing to improve people's health and fuels our community.

If you like your FITAGE product, please share your experience with us!

1. Go to Amazon / Bol.com / fitage.nl
2. Click "Your Orders" and search for the FITAGE product.
3. Click on "Write A Review" and tell us about your experience!

## ANY PROBLEMS?

If you encounter any problems, please contact us via Email.

We are available  
on Monday - Friday  
09.00-17.00 (+1 GTM)  
in English or Dutch.

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JOYTECH Healthcare Co., Ltd.
Shanghai International Holding Corp. GmbH (Europe)